



ModLink™

DICOM Conformance Statement

March 2018

LM-ENG-ML-DCS-REVB

Hyland - Proprietary

Copyright 2018 - Hyland Software, Inc.
All rights reserved.

Any comments or questions regarding the contents of this document
should be directed to the author.

Revision History

Date	Revision	Author(s)	Description
05/18/2015	A	Chris Leitner	Created.
04/2/2018	B	Ken Congdon	Final proofreading.

Table of Contents

- 1 Introduction 6**
- 2 Implementation Model..... 6**
 - 2.1 Application Data Flow Diagram 6
 - 2.2 Functional Definition of AE's 7
 - 2.3 Sequencing of Real-World Activities 7
- 3 AE Specifications 7**
 - 3.1 ModLink AE Specifications 7
 - 3.1.1 Association Establishment Policies 8
 - 3.1.2 Association Acceptance Policy 8
- 4 Communication Profiles 12**
 - 4.1 Supported Communication Stacks 12
 - 4.2 TCP/IP Stack 12
 - 4.2.1 Physical Media Support 12
- 5 Extension/Specialization/Privatization 12**
- 6 Configuration 12**
- 7 Extended Character Sets 12**

List of Figures

Figure 1. ModLink Implementation Model 6

Figure 2. Implementation Identifying Information 8

1 Introduction

This conformance statement is designed to communicate technical information regarding the ModLink product and its compliance to the DICOM 3.0 standard. ModLink provides users a simple method of receiving DICOM structured reports and forwarding that data to a voice dictation system.

2 Implementation Model

2.1 Application Data Flow Diagram




Figure 1. ModLink Implementation Model

ModLink provides users a simple method of receiving DICOM structured reports and forwarding that data to a voice dictation system.

2.2 Functional Definition of AE's

The ModLink Application Entity supports the following SCP functions:

- ▶ Receive Images This AE provides the ability to receive and store images as a DICOM Storage SCP.
- ▶ Verify Connectivity This AE provides the ability to acknowledge DICOM network connectivity as a DICOM Verification SCP.

2.3 Sequencing of Real-World Activities

Not applicable.

3 AE Specifications

3.1 ModLink AE Specifications

The ModLink AE provides standard conformance to the following DICOM 3.0 SOP classes as an SCP. Please, note ModLink does not store any received image files, image SOPs are only supported to allow for modalities that cannot send SR files only.

SOP Class Name	SOP Class UID
Verification	1.2.840.10008.1.1
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3
CR Image Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
DX Image Storage (Presentation)	1.2.840.10008.5.1.4.1.1.1.1
DX Image Storage (Raw)	1.2.840.10008.5.1.4.1.1.1.1.1
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
US Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1
US Multi-Frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3
US Multi-Frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
SC Image Storage	1.2.840.10008.5.1.4.1.1.7
MG Storage (Presentation)	1.2.840.10008.5.1.4.1.1.1.2
MG Storage (Raw)	1.2.840.10008.5.1.4.1.1.1.2.1
Multi-frame True Color Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7.4
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2
NM Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5
NM Image Storage	1.2.840.10008.5.1.4.1.1.20
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1
RF Image Storage	1.2.840.10008.5.1.4.1.1.12.2
Standard PET Image	1.2.840.10008.5.1.4.1.1.128
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33

3.1.1 Association Establishment Policies

3.1.1.1 General

The maximum PDU size for any association establishment that is offered is 512 Kbytes.

3.1.1.2 Number of Associations

The ModLink AE can establish up to twenty simultaneous associations.

3.1.1.3 Asynchronous Nature

The ModLink AE does not support asynchronous communication.

3.1.1.4 Implementation Identifying Information

The implementation identifying information for this DICOM 3.0 implementation is:

Implementation Class UID	1.3.6.1.4.1.23849.1
Version Name	PACSGEAR_v3

Figure 2. Implementation Identifying Information

3.1.2 Association Acceptance Policy

3.1.2.1 Real-World Activity – Verify Connectivity

The ModLink AE will accept associations for C-Echo and provide standard conformance to the DICOM Verification Service class.

3.1.2.1.1 Proposed Presentation Contexts

Presentation Context Table				
Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Verification	1.2.840.10008.1.2	Implicit VR Little Endian	SCP	None

3.1.2.2 Real-world Activity – Receive Images

3.1.2.2.1 The ModLink AE will accept associations for C-Storage requests and provide standard conformance to the DICOM Storage Service class for the purpose of data collection and forwarding.

3.1.2.2.2 Presentation Contexts

Proposed Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
CR Image Storage	1.2.840.10008.5.1.4.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

Proposed Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
Breast Tomo-Synthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
DX Image Storage (Presentation)	1.2.840.10008.5.1.4.1.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
DX Image Storage (Raw)	1.2.840.10008.5.1.4.1.1.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
US Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

Proposed Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
US Multi-Frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
US Multi-Frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
SC Image Storage	1.2.840.10008.5.1.4.1.1.7	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Multi-frame True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
MG Storage (Presentation)	1.2.840.10008.5.1.4.1.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
MG Storage (Raw)	1.2.840.10008.5.1.4.1.1.1.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
NM Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

Proposed Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
NM Image Storage	1.2.840.10008.5.1.4.1.1.20	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
RF Image Storage	1.2.840.10008.5.1.4.1.1.12.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Standard PET Image	1.2.840.10008.5.1.4.1.1.128	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
Comp. SR Storage	1.2.840.10008.5.1.4.1.1.88.33	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

4 Communication Profiles

4.1 Supported Communication Stacks

The ModLink AE provides DICOM 3.0 TCP/IP network communication support as defined in PS 3.8.

4.2 TCP/IP Stack

The ModLink AE implements DICOM 3.0 on top of the Windows TCP/IP stack.

4.2.1 Physical Media Support

The ModLink AE is indifferent to the physical medium over which TCP/IP executes.

5 Extension/Specialization/Privatization

Not applicable.

6 Configuration

The following items related to DICOM are configurable for the ModLink AE:

- ▶ Local AE Title
- ▶ Local Ports
- ▶ Server Socket Timeout

7 Extended Character Sets

Not applicable.