

PACSGEAR WEB FORMS



Eliminate paper from radiology reporting



Access forms on any computer, tablet or HTML 5 browser



Map SR data and technologist observations to dictation reports

Web Forms from Hyland Healthcare is a true, thin-client electronic forms solution that enables additional clinical and historical data to be mapped along with DICOM Structured Report (SR) data from PACSgear ModLink into dictation reporting systems.

The Web Forms solution allows the user to eliminate paper by creating and managing unique electronic forms that mimic their existing paper forms. Users can access the forms from anywhere on any computer or iOS, Android or Windows tablet via a URL hosted in your IT department. The user can then electronically send the form to the PACS study as an image. Furthermore, any data in the Web Form can be mapped to a custom field in the radiology report template as required.

DICOM SR data mapped to a voice dictation system enhances report productivity and quality. However, there often is other historical or clinical information not included in a DICOM SR. Technologist observations such as fetal elevation, placenta position, placental grade, murphy sign, amniotic fluid, patient history, date of injury and more are not normally part of an SR, but important clinical observations for the radiologist's interpretation. Using Web Forms, any of these observations can now be mapped to a dictation report if desired.

When integrated with ModLink, PACSgear Web Forms adds tremendous value to DICOM SR integrations by ensuring the entire scope of patient information is captured by dictation reporting systems while eliminating handwritten paper documentation. All of this adds clinical value because radiologists have complete content at their fingertips to make the most informed interpretations while enhancing productivity and workflow.

WEB FORMS HIGHLIGHTS:

- Eliminate all paper
- Integrates seamlessly with PACSgear ModLink to provide complete capture of DICOM SR and associated information
- Duplicate any paper form into a Web Forms electronic format
- Map any DICOM SR or manually entered data directly into a voice recognition/dictation system
- Launch the Web Form in context from PACS
- The Web Form can be electronically sent to the PACS System as a DICOM image for reference
- Utilize pull-down menus, check boxes, free-text fields and anatomical images
- Make key fields mandatory
- Save forms
- Create two-page forms
- Operates on any HTML 5 web browser - iOS (iPad), Android, Windows systems, etc.
- Send any DICOM SR or Web Form data to EHR or other database



Hyland[®]

With Web Forms, radiologists have complete content at their fingertips to make the most informed interpretations.

Web Forms is a simple 1-2-3 process:

1. Create a JPG, PNG or GIF image of your paper form and upload to Web Forms
2. Open the form image in Web Forms to use as a background
3. Insert data fields anywhere you want from DICOM Modality Worklist (i.e. arrow at top with Patient Name, DOB, etc.), measurement fields (arrow at bottom) from DICOM SR, or manually enter data.

Learn more at HylandHealthcare.com/EnterpriseImaging

Hyland ModLink Web Forms
 Demo OB Single

Hyland HEALTHCARE Name: _____
 DOB _____
 Accession # _____
 MR # _____
 Date: _____

Sonographer worksheet only:
 not for diagnostic purposes

Single Complete OB Ultrasound

Symptoms _____
 Sonographer ADR Phone _____

Comparison <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	1st US _____ Most Recent US _____
Transvaginal scan <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	
Fetal Presentation <input type="checkbox"/>	
Fetal Movement <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	
Cardiac Rate <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	rate _____
Amniotic Fluid <input checked="" type="checkbox"/> Normal <input type="checkbox"/> Abnormal	
Placenta Position <input type="checkbox"/>	

LMP <input type="checkbox"/> W <input type="checkbox"/> D <input type="checkbox"/> EDD	<input type="checkbox"/> Unknown
GA by 1st US <input type="checkbox"/> W <input type="checkbox"/> D <input type="checkbox"/> EDD	Should be today
GA by today's US <input type="checkbox"/> W <input type="checkbox"/> D <input type="checkbox"/> EDD	
During by <input type="checkbox"/> LMP <input type="checkbox"/> 1st US <input type="checkbox"/> IVF/ARI	

Comments/Other Findings

not a chart form

			EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands	Australian Sponsor Emergo Australia Level 20 Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia
		Hyland LLC 4309 Hacienda Drive, Suite 500 Pleasanton, California 94588 USA	LX-PMK-MLF-BRO-EN-REVC	

