

PACSGEAR MEDIA WRITER

- ✓ Easy to configure, deploy and maintain
- ✓ Burn DICOM images, PDF and reports
- ✓ Intuitive client software lets you work anywhere

Affordable, reliable DICOM CD/DVD burning

Media Writer from Hyland Healthcare delivers an affordable, reliable DICOM CD/DVD/USB burning solution for images, regardless of origin. Media Writer provides a media distribution system that is easy to configure, deploy and maintain for imaging-based studies such as radiology, cardiology, gastroenterology, dermatology, orthopedics, digital pathology and dental.


Media Writer is a low-maintenance, cost-effective and intuitive software/hardware solution that allows your organization to provide this longstanding, important distribution method for providers and patients alike.

MEDIA WRITER FEATURES


- Directly burns or uploads DICOM CDs/DVDs from any VNA/PACS archive or modality user interface; intuitive client software lets you work anywhere.
- Safeguards your portable media with AES 256-bit encryption technology.
- For large studies, Media Writer spans multiple CDs or autodetects between CDs and DVDs.
- Includes Gear View Basic, an embedded DICOM viewer.
- Exports and anonymizes DICOM CDs/DVDs with one click.
- Creates customized splash screens and CD/DVD labels.
- Combines studies from multiple DICOM devices.
- Supports multipatient burning (i.e. conferences, Tumor Board).
- Interface supports 12 languages.
- Supports IHE workflow.

The Hyland logo is positioned in the bottom left corner of the slide. It consists of the word "Hyland" in a white, serif font, centered within a square that has a green-to-blue gradient background.

MEDIA WRITER TECHNICAL SPECIFICATIONS



Model	Media WriterD200	Media WriterD100
Burner/printer	Epson® Discproducer™	Epson® Discproducer™
Robotics	Yes, AcuGrip™ technology	Yes, AcuGrip™ technology
Typical usage	Hospitals, film rooms (up to 200 discs per day)	Departments, imaging centers, medical office buildings (up to 100 discs per day)
Burn/print time per disc	3-4 minutes, typical	3-4 minutes, typical
Disc capacity	Holds 100 discs (CD and/or DVD)	Holds 100 discs (CD and/or DVD)
CD/DVD drives	2 CD/DVD drives	1 CD/DVD drive
Recording speed	40x CD-R; 12x DVD3R; 8x DVD3R DL	40x CD-R; 12x DVD3R; 8x DVD3R DL
Media	Inkjet-printable 12cm	Inkjet-printable 12cm
Inkjet print method	MicroPiezo® technology	MicroPiezo® technology
Ink cartridge/type	Dye Ink/6 cartridge	Dye Ink/6 cartridge
Ink colors	Cyan, magenta, yellow, light cyan, light magenta, black	Cyan, magenta, yellow, light cyan, light magenta, black
Palette yield	16.7M colors, 1000 discs per six cartridges, typical	16.7M colors, 1000 discs per six cartridges, typical
Data interface	USB 3.0 high speed	USB 3.0 high speed
OS support	Windows 10 x32/64	Windows 10 x32/64
Electrical	100-120VAC/240VAC, 50/60Hz, 50W	100-120VAC/240VAC, 50/60Hz, 50W
Certifications	FCC Class B, CAN/CSA Class B, CE	FCC Class B, CAN/CSA Class B, CE
System requirements	1.4GHz processor, 10GB hard drive space, 1GB memory	1.4GHz processor, 10GB hard drive space, 1GB memory
Dimensions (HxWxD)	14" x 15" x 18" (348 x 377 x 465mm)	12.76" x 14.84" x 18.31" (324 x 377 x 465mm)
Weight	52.9 lbs. (24kg)	46.2 lbs. (21kg)


 EU Authorized Representative
EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands


Manufacturer
Hyland LLC
 4309 Hacienda Drive, Suite 500
 Pleasanton, CA 94588 USA

For support, contact:

Tel: 1-844-535-1404

Web: <https://community.hyland.com>

Australian Sponsor:
 Emergo Australia
 201 Sussex Street, Level 20,
 Tower II Darling Park,
 Sydney, NSW 2000, Australia

Learn more at HylandHealthcare.com/EnterpriseImaging

Hyland™

Media Writer is not used for the purpose of diagnosis, cure, mitigation, prevention, or treatment of a disease or condition.

Part number: DOC-MW-BRO-EN-Rev A

©2022 Hyland Software, Inc. and its affiliates. All rights reserved. All Hyland product names are registered or unregistered trademarks of Hyland Software, Inc. or its affiliates.