

PACS SCAN

✓ Easily scan documents and create electronic forms from any department to any PACS/EMR

✓ Quickly import and localize images (JPEG/AVI/MPEG) and DICOM CDs/DVDs

✓ Compatible with more than 25 Hyland solutions and TWAIN compliant scanners

Industry-leading document and multimedia connectivity software

PACS Scan from Hyland Healthcare lets you scan documents, create electronic forms, import JPEGs, AVIs and MPEGs from any department to any picture and archiving communication system (PACS)/electronic medical record (EMR). Convert and send Word files, PDFs, reports, emails and other images with the PACS Scan Virtual Printer, a printer selection that appears within your Windows application.

Eliminate paper, improve workflow and provide access to important clinical documentation, such as tech notes, screening forms and order information. In PACS/EMR, view documents, electronic forms and clinical images along with patient studies at the point of care.

PACS Scan includes additional features such as electronic forms and MediaImport to import and localize DICOM CDs/DVDs from a wide range of PACS vendors.

PACS SCAN FEATURES

- Simple user interface allows anyone to scan documents without prior knowledge of PACS. No training required.
- PACS Scan Virtual Printer is now integrated with PACS Scan.
- Interfaces to radiological information system (RIS) worklist or DICOM Archive to provide automatic loading of patient demographics for ease of use.
- DICOM Modality Worklist (MWL) ensures an accurate match with patient studies.
- Includes extensive bar code support for rapid demographic entry.
- Scan and send documents to PACS/EMR using TWAIN-compliant scanners.
- Validated interoperability with over 25 PACS.
- Direct import of multimedia formats including JPEG, bitmap, TIFF, PNG, AVI/MPEG and DICOM CD/DVD.
- Drag and drop functionality is available to import JPEG/AVI/MPEG images.
- MediaImport (DICOM CD/DVD import) supports automatic decompression of compressed JPEG, JPEG2000 or RLE DICOM studies before sending them.
- Access to confidential patient information is protected by username and password for each user.
- Supports LDAP/Active Directory for single sign-on convenience.




Hyland[®]

PACS SCAN

1. Select patient/study demographics from your RIS or DICOM archive.




2. Scan the document and click "send."




3. The document is now matched with the appropriate study.

PACS SCAN VIRTUAL PRINTER



4. From your Windows-based application, select File > Print, and then select "PacsScan" as the printer.

5. Select patient/study demographics from your RIS or DICOM archive, click "Preview" to preview the document or click "Send" to convert and send the document to PACS/EMR.



6. The document is now matched with the appropriate study.



Rapid data entry with stylus tablet

ELECTRONIC FORMS

PACS Scan allows you to create and distribute electronic forms to improve efficiency and workflow, which eliminates paper scanning entirely!

- Eliminate paper "Print & Scan" cycle.
- Use optional stylus tablet for rapid data entry.
- Auto-fill key fields with MWL.

DICOM CD/DVD IMPORT

- Reads all DICOM-compliant and some not-so-DICOM-compliant CD formats.
- Imports patient demographics — supports DICOM Modality Worklist for accurate demographic entry.
- Edits and localizes patient demographics — changes fields such as Patient ID/MRN and accession number, prior to import.

EU Authorized Representative
EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands

Manufacturer
Hyland LLC
 4309 Hacienda Drive, Suite 500
 Pleasanton, CA 94588 USA

Australian Sponsor:
 Emergo Australia
 201 Sussex Street, Level 20,
 Tower II Darling Park,
 Sydney, NSW 2000, Australia

For support, contact:
 Tel: 1-844-535-1404
 Web: <https://community.hyland.com>

Learn more at [HylandHealthcare.com](https://www.hylandhealthcare.com)

PACS Scan is not used for the purpose of diagnosis, cure, mitigation, prevention, or treatment of a disease or condition.

Part number: DOC-PS-BRO-EN-RevA

©2020 Hyland Software, Inc. and its affiliates. All rights reserved. Hyland Healthcare is part of Hyland Software, Inc. and its affiliated companies. All Hyland product names are registered or unregistered trademarks of Hyland Software, Inc. or its affiliates in the United States and other countries.

Hyland